

Universal FFPE IHC Tissue Microarray

GENERAL INFORMATION

Product Name: Universal FFPE IHC Tissue Microarray

Reference Number: 7310-9120 5-pack Slides

Date of Manufacturing: See Product Label

Intended Use: For Research Use Only

DESCRIPTION

Tissue info: Tissue See Array Map

of Origin: Organ Human

type: Fixation Breast, Tonsil, Appendix, Placenta

Condition: 10% neutral buffered formalin (NBF) for 24 hours at 24-27°C

Tissue Size: Approximately 2mm wide, 10 mm long and greater than 1mm thick.

STORAGE AND STABILITY

Storage Condition	Stability*
Ambient temperature	6 months
2-8 ° C with desiccation (Recommended)	2 years
-20 ° C to -10 ° C	5 years

*Guideline for general applications, such as immunohistochemistry (IHC) or DNA in situ hybridization (ISH). Certain biomolecules may be less stable during storage.

SAFETY AND PRECAUTIONS

This product does not contain hazardous material. Wear appropriate personal protective equipment (PPE) when handling reagents and biological specimens.

RECOMMENDED PROCEDURES


Staining using FFPE slides:

1. Bake slides at 60°C for 30-60 min.
2. Deparaffinize two times in Xylene or Xylene substitute for 5 min each time.
3. Rinse two times with 100% ethanol for 1 min each time.
4. Air dry slides or rehydrate in ethanol series (95% 2 min, 70%, 2 min, 50% 2 min, 1XPBS 2 min).
5. Proceed to staining protocol.



Acepix Biosciences, Inc.
32990 Alvarado-Niles Road Ste 910
Union City · CA · 94587 · USA
Tel: 1 (888) 866-0580
Web: www.acepixon.com

ARRAY MAP



	1	2	3	4
A	Hs Normal Breast	Hs Normal Tonsil	Hs Normal Appendix	Hs Normal Placenta