

A-FLX™ Animal Tissue FFPE BALB/C Mouse Spleen

GENERAL INFORMATION

| | |
|------------------------|--|
| Product Name: | A-FLX™ Animal Tissue FFPE BALB/C Mouse Spleen |
| Reference Number: | 7010-1410 Block 7010-1420 Slide (5µm) |
| Date of Manufacturing: | See product label |
| Lot Number: | See product label |
| Intended Use: | For research use only |

DESCRIPTION

| | |
|---------------------|--|
| Species: | Mouse |
| Strain: | BALB/C |
| Tissue: | Spleen |
| Fixation Condition: | 10% neutral buffered formalin (NBF) for 24 hours at 24-27°C |
| Product Format: | |
| Block: | Paraffin embedded block |
| Slide: | One unstained section mounted on Superfrost™ Plus slide. Section thickness: 5µm |

STORAGE AND STABILITY

| Storage Condition | Stability* |
|---------------------------------------|------------|
| Ambient temperature | 6 months |
| 2-8 °C with desiccation (Recommended) | 2 years |
| -20 °C to -10 °C | 5 years |

*Guideline for general applications, such as immunohistochemistry (IHC) or DNA *in situ* hybridization (ISH). Certain biomolecules may be less stable during storage.

SAFETY AND PRECAUTIONS

This product does not contain hazardous material. Wear appropriate personal protective equipment (PPE) when handling reagents and biological specimens.



Acepix Biosciences, Inc.

31351 Medallion Drive · Hayward · CA · 94544 · USA
Tel: 1 (888) 866-0580 · Web: www.acepixon.com

