

A-FLX™ Animal Tissue FFPE BALB/C Mouse Bladder

GENERAL INFORMATION

Product Name:	A-FLX™ Animal Tissue FFPE BALB/C Mouse Bladder
Reference Number:	7010-2210 Block 7010-2220 Slide (5µm)
Date of Manufacturing:	See product label
Lot Number:	See product label
Intended Use:	For research use only

DESCRIPTION

Species:	Mouse
Strain:	BALB/C
Tissue:	Bladder
Fixation Condition:	10% neutral buffered formalin (NBF) for 24 hours at 24-27°C
Product Format:	
Block:	Paraffin embedded block
Slide:	One unstained section mounted on Superfrost™ Plus slide. Section thickness: 5µm

STORAGE AND STABILITY

Storage Condition	Stability*
Ambient temperature	6 months
2-8 °C with desiccation (Recommended)	2 years
-20 °C to -10 °C	5 years

*Guideline for general applications, such as immunohistochemistry (IHC) or DNA *in situ* hybridization (ISH). Certain biomolecules may be less stable during storage.

SAFETY AND PRECAUTIONS

This product does not contain hazardous material. Wear appropriate personal protective equipment (PPE) when handling reagents and biological specimens.



Acepix Biosciences, Inc.

31351 Medallion Drive · Hayward · CA · 94544 · USA
Tel: 1 (888) 866-0580 · Web: www.acepixon.com