

## A-FLX™ Animal Tissue FFPE BALB/C Mouse Adrenal Gland

### GENERAL INFORMATION

Product Name: A-FLX™ Animal Tissue  
FFPE BALB/C Mouse Adrenal Gland

Reference Number: 7010-2410 Block  
7010-2420 Slide (5µm)

Date of Manufacturing: See product label

Lot Number: See product label

Intended Use: For research use only

### DESCRIPTION

Species: Mouse

Strain: BALB/C

Tissue: Adrenal Gland

Fixation Condition: 10% neutral buffered formalin (NBF) for 24 hours at 24-27°C

Product Format:

Block: Paraffin embedded block

Slide: One unstained section mounted on Superfrost™ Plus slide. Section thickness: 5µm

### STORAGE AND STABILITY

Storage Condition	Stability*
Ambient temperature	6 months
2-8 °C with desiccation (Recommended)	2 years
-20 °C to -10 °C	5 years

\*Guideline for general applications, such as immunohistochemistry (IHC) or DNA *in situ* hybridization (ISH). Certain biomolecules may be less stable during storage.

### SAFETY AND PRECAUTIONS

This product does not contain hazardous material. Wear appropriate personal protective equipment (PPE) when handling reagents and biological specimens.



#### Acepix Biosciences, Inc.

31351 Medallion Drive · Hayward · CA · 94544 · USA  
Tel: 1 (888) 866-0580 · Web: [www.acepixon.com](http://www.acepixon.com)

