

A-FLX™ Animal Tissue FFPE C57BL/6 Mouse Kidney

GENERAL INFORMATION

Product Name:	A-FLX™ Animal Tissue FFPE C57BL/6 Mouse Kidney
Reference Number:	7011-0810 Block 7011-0820 Slide (5µm)
Date of Manufacturing:	See product label
Lot Number:	See product label
Intended Use:	For research use only

DESCRIPTION

Species:	Mouse
Strain:	C57BL/6
Tissue:	Kidney
Fixation Condition:	10% neutral buffered formalin (NBF) for 24 hours at 24-27°C
Product Format:	
Block:	Paraffin embedded block
Slide:	One unstained section mounted on Superfrost™ Plus slide. Section thickness: 5µm

STORAGE AND STABILITY

Storage Condition	Stability*
Ambient temperature	6 months
2-8 °C with desiccation (Recommended)	2 years
-20 °C to -10 °C	5 years

*Guideline for general applications, such as immunohistochemistry (IHC) or DNA *in situ* hybridization (ISH). Certain biomolecules may be less stable during storage.

SAFETY AND PRECAUTIONS

This product does not contain hazardous material. Wear appropriate personal protective equipment (PPE) when handling reagents and biological specimens.



Acepix Biosciences, Inc.

31351 Medallion Drive · Hayward · CA · 94544 · USA
Tel: 1 (888) 866-0580 · Web: www.acepixon.com