

## A-FLX™ Animal Tissue FFPE C57BL/6 Mouse Pancreas

### GENERAL INFORMATION

Product Name:	A-FLX™ Animal Tissue FFPE C57BL/6 Mouse Pancreas
Reference Number:	7011-1210 Block 7011-1220 Slide (5µm)
Date of Manufacturing:	See product label
Lot Number:	See product label
Intended Use:	For research use only

### DESCRIPTION

Species:	Mouse
Strain:	C57BL/6
Tissue:	Pancreas
Fixation Condition:	10% neutral buffered formalin (NBF) for 24 hours at 24-27°C
Product Format:	
Block:	Paraffin embedded block
Slide:	One unstained section mounted on Superfrost™ Plus slide. Section thickness: 5µm

### STORAGE AND STABILITY

Storage Condition	Stability*
Ambient temperature	6 months
2-8 °C with desiccation (Recommended)	2 years
-20 °C to -10 °C	5 years

\*Guideline for general applications, such as immunohistochemistry (IHC) or DNA *in situ* hybridization (ISH). Certain biomolecules may be less stable during storage.

### SAFETY AND PRECAUTIONS

This product does not contain hazardous material. Wear appropriate personal protective equipment (PPE) when handling reagents and biological specimens.



#### Acepix Biosciences, Inc.

31351 Medallion Drive · Hayward · CA · 94544 · USA  
Tel: 1 (888) 866-0580 · Web: [www.acepixon.com](http://www.acepixon.com)