

## A-FLX™ Animal Tissue FFPE C57BL/6 Mouse Eye

### GENERAL INFORMATION

Product Name: A-FLX™ Animal Tissue  
FFPE C57BL/6 Mouse Eye  
Reference Number: 7011-1510 Block  
7011-1520 Slide (5µm)  
Date of Manufacturing: See product label  
Lot Number: See product label  
Intended Use: For research use only

### DESCRIPTION

Species: Mouse  
Strain: C57BL/6  
Tissue: Eye  
Fixation Condition: 10% neutral buffered formalin (NBF) for 24 hours  
at 24-27°C  
Product Format:  
Block: Paraffin embedded block  
Slide: One unstained section mounted on Superfrost™  
Plus slide. Section thickness: 5µm

### STORAGE AND STABILITY

Storage Condition	Stability*
Ambient temperature	6 months
2-8 °C with desiccation (Recommended)	2 years
-20 °C to -10 °C	5 years


\*Guideline for general applications, such as immunohistochemistry (IHC) or DNA *in situ* hybridization (ISH). Certain biomolecules may be less stable during storage.

### SAFETY AND PRECAUTIONS

This product does not contain hazardous material. Wear appropriate personal protective equipment (PPE) when handling reagents and biological specimens.



#### Acepix Biosciences, Inc.

 31351 Medallion Drive · Hayward · CA · 94544 · USA  
Tel: 1 (888) 866-0580 · Web: [www.acepixon.com](http://www.acepixon.com)