

A-FLX™ Animal Tissue FFPE C57BL/6 Mouse Intestine

GENERAL INFORMATION

Product Name: A-FLX™ Animal Tissue
FFPE C57BL/6 Mouse Intestine

Reference Number: 7011-1810 Block
7011-1820 Slide (5µm)

Date of Manufacturing: See product label

Lot Number: See product label

Intended Use: For research use only

DESCRIPTION

Species: Mouse

Strain: C57BL/6

Tissue: Intestine

Fixation Condition: 10% neutral buffered formalin (NBF) for 24 hours at 24-27°C

Product Format:

Block: Paraffin embedded block

Slide: One unstained section mounted on Superfrost™ Plus slide. Section thickness: 5µm

STORAGE AND STABILITY

Storage Condition	Stability*
Ambient temperature	6 months
2-8 °C with desiccation (Recommended)	2 years
-20 °C to -10 °C	5 years

*Guideline for general applications, such as immunohistochemistry (IHC) or DNA *in situ* hybridization (ISH). Certain biomolecules may be less stable during storage.

SAFETY AND PRECAUTIONS

This product does not contain hazardous material. Wear appropriate personal protective equipment (PPE) when handling reagents and biological specimens.



Acepix Biosciences, Inc.

31351 Medallion Drive · Hayward · CA · 94544 · USA
Tel: 1 (888) 866-0580 · Web: www.acepixon.com