

A-FLX™ Animal Tissue FFPE Wistar Rat Eye

GENERAL INFORMATION

Product Name: A-FLX™ Animal Tissue
FFPE Wistar Rat Eye
Reference Number: 7020-1510 Block
7020-1520 Slide (5µm)
Date of Manufacturing: See product label
Lot Number: See product label
Intended Use: For research use only

DESCRIPTION

Species: Rat
Strain: Wistar
Tissue: Eye
Fixation Condition: 10% neutral buffered formalin (NBF) for 24 hours at 24-27°C
Product Format:
Block: Paraffin embedded block
Slide: One unstained section mounted on Superfrost™ Plus slide. Section thickness: 5µm

STORAGE AND STABILITY

| Storage Condition | Stability* |
|---------------------------------------|------------|
| Ambient temperature | 6 months |
| 2-8 °C with desiccation (Recommended) | 2 years |
| -20 °C to -10 °C | 5 years |

*Guideline for general applications, such as immunohistochemistry (IHC) or DNA *in situ* hybridization (ISH). Certain biomolecules may be less stable during storage.

SAFETY AND PRECAUTIONS

This product does not contain hazardous material. Wear appropriate personal protective equipment (PPE) when handling reagents and biological specimens.



Acepix Biosciences, Inc.

31351 Medallion Drive · Hayward · CA · 94544 · USA
Tel: 1 (888) 866-0580 · Web: www.acepixon.com

